



Inspira – Advanced Safety Ventilator (ASV)

Breakthrough

New



- Ventilates mice to cats with one ventilator (15 g to 10 kg)
- Easy to use
- Two models available:
 - Volume-controlled
 - Pressure-controlled
- Airway pressure monitoring detects over/under pressurization of lungs
- Assist Mode helps wean animal off the ventilator
- Built in sigh breath maintains optimal lung functioning
- Adjustable I:E ratio (1:4 to 4:1)
- Ability to ventilate with air or non-flammable

Newer Information
About This Product
Is Located
On Page 174
In Our
New Products Catalog
[Click Here](#)

Announcing two new ad

- Inspira Advanced S
- Inspira Advanced S

Both ventilators are the r
ment. Each features an e
display, and all gas ports
functions include Airway
Programmable Sigh Bre

Coming soon, the Inspira
The ARV builds on the ASV
of optional features inclu

- Sidestream O₂ mo
- Mainstream or si
- Sidestream anesth
- Capnography
- Pulse oximetry
- Remote airway pressure monitoring

Call Harvard Apparatus for more details.

mals

ilator that can ventilate animals
chine. Respiration rates range
e adjusts from 0.1 ml to 100 ml.
e one ventilator for mice and an-
can do it all.

press run! A microprocessor rec-
a user-entered animal weight. The
meters if needed. The ventilators
d simple to use keypad.

the front of the machine for conve-
vo analog outputs, an external trig-
are also located on the front of the
n the back of the machine and can
ing, or printer connection.

g

the animal's weight and press run.
sets a standard set of respiratory
ers is quick and easy. The bright,
umeric keypad and intuitive func-

tion keys simplify parameter entry. Changing a value (e.g., respiratory rate) requires only a few keystrokes. The user is prompted at each step of the process to assure you of quick accurate setup. As a final check, each function key (rate, volume, I:E ratio, sigh and pressure) can be pressed to display the current value being used.

Safe Range™

Setting up the ASV ventilators is both easy and safe thanks to the Safe Range™ software. Once the animal's body weight has been entered the software computes the median settings for tidal volume and respiratory rate. These same equations are used to calculate a Safe Range™ around the median values. All the default parameter settings can be overridden. If you enter a setting outside of the Safe Range for the body mass of the animal to be ventilated you are warned that the value is outside the Safe Range™ and asked to confirm the setting.

